

INFUSEIT™ PRESSURE INFUSOR

LATEX FREE | DISPOSABLE AND REUSABLE | DEHP FREE | MR CONDITIONAL

InfuseIT™ pressure infusors are designed for ease of use. Faster inflation, individual hooks for solution and IV poles, and highly visible pressure gauge facilitate fast setup and easy maintenance. The enclosed valve eliminates the risk of fluids being trapped and hidden and the MRI conditional spring makes it safe to use within an MRI environment. Available in disposable and extended use reusable versions.



Clear back panel provides full visibility of solution bag to monitor fluid level and easily scan patient barcode.



Palm-fitted, efficient hand bulb makes inflation easy.



Stopcock Valve

Highly visible pressure gauge provides easy to read 360° window which can be read from a distance



IV fluid bag hook is durable and designed to fit all standard bags snugly

Reusable bags are navy color and disposable bags are white.



0-150 mm Hg



150-250 mm Hg



250-300 mm Hg



300+ mm Hg

MRI conditional spring

INFUSEIT™ PRESSURE INFUSOR

- Ergonomically designed hand bulb for fast inflation
- Pressure relief valve prevents over inflation (375 mm Hg pressure relief)
- Enclosed valve helps eliminate the risk of fluids being trapped and hidden
- Pressure gauge utilizes a non-ferrous bronze spring making it MRI conditional
- Reusable InfuseIT models constructed of nylon for extended use
- Incredibly reliable - 100% tested

DISPOSABLE INFUSEIT™ PRESSURE INFUSOR

| ITEM | DESCRIPTION | VOLUME | PK |
|-------------|-------------|---------|----|
| ZIT-500-10 | InfuseIT | 500 mL | 10 |
| ZIT-1000-10 | InfuseIT | 1000 mL | 10 |
| ZIT-3000-10 | InfuseIT | 3000 mL | 10 |

REUSABLE INFUSEIT™ PRESSURE INFUSOR

| ITEM | DESCRIPTION | VOLUME | PK |
|------------|-------------|---------|----|
| ZIT-520-5 | InfuseIT | 500 mL | 5 |
| ZIT-1020-5 | InfuseIT | 1000 mL | 5 |